



BIOTECH **SUPPORT GROUP**

NRicher™ Mx

General Enrichment for All Biofluids and Tissue Lysates


- Consumable chemically derived beads, species agnostic as they are not derived from antibodies
- Enrich low abundance proteomes from any source, from sera/plasma to cell lysates from both animals and humans, >90% Albumin removal
- Scalable protocol from small to large sample volumes, from 10 to 500 µl, and low to high protein concentrations
- Enriched sub-proteome, for better signal quantitation between samples
- Does not require any specialized instruments, just a standard microfuge
- Bead format suitable for automation compatibility, please inquire
- On-Bead digestion for LC-MS analysis, or optional elution for any functional, enzymatic, or immunoassay analysis

NRicher™ Mx employs the use of a bead cocktail, which allows for one, rather than multiple LC-MS analyses to establish dynamic range compression. **NRicher™ Mx** is thus an all-purpose proteomic enrichment product that can be used for any sample type, from biofluids to tissue lysates. It is compatible with up to 1% non-ionic detergent concentrations.

It is particularly useful for membrane proteins. Targets of over 50% of all therapeutic drugs, membrane proteins perform a variety of functions including:

- Receptors which relay information between internal and external environments
- Transport of molecules and ions across the cell membrane
- Enzymatic, and
- Adhesion

NRicher™ Mx beads provide excellent 2-3X enrichment of membrane proteins, most of which are not observable in neat serum.



Another example of **NRicher™ Mx** bead enrichment is α -Synuclein, a biomarker for Parkinsons Disease, observed by **NRicher™ Mx**, but not observed in the neat serum.

	Uniprot	NRicher™ Mx	Neat serum
SNCA	P37840	1.1E+06 total signal	Not detected



B|IOTECH SUPPORT GROUP


NRicher™ workflow

High abundance (i.e., Albumin) proteins selectively passes or voids through



the beads, concentrating and enriching sub-proteomes on the beads


Enriched sub-proteomes remain on the bead and digested using Bead-Assisted Sample Prep (BASP); protocols provided with the NRicher™ products



OR

Digest Options

Eluate from NRicher™ can be digested by any preferred method, In-Gel, In-Solution, FASP, etc.



LC-MS

Eluates from NRicher™ beads can be applied to other common analyses:

- Enzymatic/Functional assays
- 2DE
- ELISA/immunoassay

The NRicher™ Workflow. All NRicher™ beads are processed the same, using buffers and spin-filters provided with the kits. The beads are supplied as a dry powder, weighed and dispensed into the top of a spin-filter, and follows a bind/wash protocol using a standard microfuge to separate the buffer solutions from the beads. Once the NRicher™-derived sub-proteome (different for each application) is bound to the beads, a variety of options are available to the user including:

>Bead-Assisted Sample Prep (BASP™), whereby reduction, alkylation and digestion are performed on the bead-bound proteome, without the use of detergents, seamlessly integrating to LC-MS analysis,

OR

>Optional Elution to off-bead digestion (i.e., FASP), or other common functional or immunoassay analyses



BIOTECH SUPPORT GROUP

Product	Size	Total serum/plasma samples processed	Item No.
NRicher™ Mx	10 Preps	10 x (2-4) mg total protein samples	NIMX-10
NRicher™ Mx	50 Preps	50 x (2-4) mg total protein samples	NIMX-50

Processes 25-50 µl serum per prep, or 2-4 mg total protein from cell lysates (diluted to ≤0.1% non-ionic detergent). If the cell lysates are dilute, follow the cell lysate protocol, step # 2 in protocol to follow.

For best results – the lysate should be clear and free of colloidal material. We recommend first filtering through a 0.45 µm syringe, or microfuge-type filter before beginning the prep.

The centrifugation time may vary, adjust as necessary to get complete filtration through the beads.

The protocol can be scaled up or down proportionally to adjust for different volumes. The bead amount can be adjusted to accommodate more or less proteome capture.

In bold are the **NRicher™** kit components.

Items Required	10 Prep	50 Prep	Reagent
NRicher™ Mx Beads	0.25 gram	1.25 gram	Supplied
Binding Buffer NRBB (0.05M HEPES, pH 6.0)	60 ml	300 ml	Supplied
Wash Buffer NRWB (0.05M HEPES, pH 7.0)	12 ml	60 ml	Supplied
Elution Buffer NREB (0.25M Tris + 0.5M NaCl, pH 9-10)	3 ml	15 ml	Supplied
Spin-filter & tube assemblies*	10	50	Supplied
DTT, Iodoacetamide, Trypsin and Formic Acid, 50% Acetonitrile (ACN)			Not Supplied

***Additional Spin-Filters (low protein binding, 0.45 µm filter element) can be purchased separately, please inquire.**

If there are any questions about compatibility or substitution with other buffers, please contact us.

Protocol For Enrichment of Low Abundance Proteins, & On bead Digestion For LC-MS Analysis

Optional Elution Protocol is included for Off-bead digestion or any functional, enzymatic, or immunoassay analysis



BIOTECH SUPPORT GROUP

1. SAMPLE PROCESSING.

FOR SERUM OR BIOFLUIDS.

1a. Weigh out 25 mg of **NRicher™ Mx** beads in a spin-filter. Add 150 µl of **Binding Buffer NRBB**. Vortex for 5 minutes at room temperature followed by centrifugation for 2 minutes at 5,000-10,000 rpm (2,000-8,000xg). Discard the filtrate. Repeat step-1.

1b. Add 200 µl of **Binding Buffer NRBB** to **NRicher™** beads followed by (25 to 50) µl Serum, Plasma or other biofluid, to the beads. Vortex or mix thoroughly for 10 min and then centrifuge for 4 minutes at 5,000-10,000 rpm (2,000-8,000xg). Discard the filtrate.

FOR CELL LYSATES UP TO 5 ML (≤0.1% NON-IONIC DETERGENT), CONTAINING 2-4 MG TOTAL PROTEIN.

1a. To clarified cell lysate, add 5 ml **Binding Buffer NRBB** in 15 ml centrifuge tube.

1b. Add 25 mg **NRicher™** beads and vortex for 25 minutes. Note: Vortex sufficiently so that the beads do not settle at the bottom of the centrifuge tube.

1c. Allow the beads to settle for 10 minutes. Decant or pipette off the supernatant.

1d. Using a wide bore pipette, transfer **NRicher™** beads to the supplied spin-filter tube assembly. Note: if all beads do not transfer, use additional **Binding Buffer NRBB** (approximately 1 ml) to resuspend & transfer again.

1e. Centrifuge for 4 minutes at 5,000-10,000 rpm (2,000-8,000xg). Discard the filtrate.

2. To the **NRicher™** beads, add 250 µl of **Wash Buffer NRWB**. Vortex for 5 min and centrifuge for 4 minutes at 5,000-10,000 rpm (2,000-8,000xg). Discard the **Wash** filtrate.
3. Repeat Wash Step-2.
4. **After discarding the wash from step 4, the NRicher™ beads contain the enriched sub-proteome. As an option for LC-MS sample preparation, the bead assisted on-bead digestion protocol (BASP™) is provided starting on step 6, see box below.**

OPTIONAL BEAD ELUTION. To the beads, add 300 µl of **Elution Buffer NREB**. Vortex or mix thoroughly for 10 min and centrifuge for 4 minutes at 5,000-10,000 rpm (2,000-8,000xg). Recover the filtrate as the eluted sub-proteome (0.25M Tris + 0.5M NaCl, pH 9.0-10.0), suitable for further analysis.



B:OTECH SUPPORT GROUP

The bead assisted on-bead digestion protocol (BASP™) is provided below. The digest buffer is **Wash Buffer NRWB** (0.05M HEPES, pH 7.0). Comparable buffers (0.02-0.10M, pH 6-7) can be used. Higher pH buffers are not recommended.

5. Using **Wash Buffer NRWB**, prepare to 10mM of DTT concentration, and add 100 μ l to the **NRicher™** beads and vortex for 10 minutes and incubate for 30 minutes at 60C.
6. Cool the samples to RT, add suitable volume of Iodoacetamide to 20mM and incubate in the dark for 45 minutes.
7. Centrifuge at 5,000-10,000 rpm (2,000-8,000xg) for 5 minutes, and discard filtrate. Rinse the bottoms of the spin-filter tubes with 500 μ l of 50% ACN, **Wash Buffer NRWB** twice, to remove any traces of the filtrate.
8. Add 8 μ g trypsin in 100 μ l **Wash Buffer NRWB** to the **NRicher™** beads and keep at 37°C for a minimum 4 hours to maximum overnight. Overnight is recommended to start with. In select targeted circumstances, 2 hours may be sufficient.
9. Centrifuge at 5,000-10,000 rpm (2,000-8,000xg) for 5 minutes, and retain digested peptides filtrate.
10. To further extract remaining peptides, add 150 μ L 10% formic acid, vortex 10 min, centrifuge at 5,000-10,000 rpm (2,000-8,000xg) for 5 mins., and combine this volume with volume from step 10.
11. Total is about 250 μ l. Prepare to desired final concentration. Store at -80°C until LC-MS/MS.



BIOTECH SUPPORT GROUP

NRicher™ Beads Are Versatile to A Variety of Bead Processing Formats

Microfuge Spin-filter is our standard



Other formats compatible with the 50 µm NRicher™ beads are:

High Throughput Automation Compatible INTip™ SPE (DPX Technologies) Format

Aspirate and dispense
cycles mix NRicher™ beads
and solutions



The INTip™ SPE tip format improves ease of use and scalability to process multiple samples in parallel, utilizing 96-well plates and automated liquid handlers. The tip-based formats have been proven to be compatible with most automation platforms, i.e., Integra, Hamilton, etc. Please inquire for more information, as these formats are customized to the application and automation platform.

96-Well Vacuum or Pressure Filter Format

The NRicher™ beads can be readily processed in 96-well filter formats. Please inquire.





BIOTECH SUPPORT GROUP

Related Separations, Enrichment/Depletion & Sample Prep - All Product Categories

(https://www.biotechsupportgroup.com/Products-a-z_a/258.htm)

Albumin & IgG Removal

(<https://www.biotechsupportgroup.com/Articles.asp?ID=451>)

Lipid Removal and Clarification

(<https://www.biotechsupportgroup.com/Articles.asp?ID=456>)

Hemoglobin Removal

(<https://www.biotechsupportgroup.com/Articles.asp?ID=452>)

Sample Prep – Liquid Biopsy

(<https://www.biotechsupportgroup.com/Articles.asp?ID=457>)

Sample Prep – Glyco, Virus, Kinase, Aqueous Protein Crash/Metabolomics

(<https://www.biotechsupportgroup.com/Articles.asp?ID=453>)

Sample Prep – Mass Spectrometry

(<https://www.biotechsupportgroup.com/Articles.asp?ID=432>)

Sample Prep – Genomics

(<https://www.biotechsupportgroup.com/Articles.asp?ID=455>)

CONTACT US

We welcome your questions and comments regarding our products.

Tel: 732-274-2866, 800-935-0628 (North America) Mon – Fri 9am-6pm EST.

Email: sales@biotechsupportgroup.com

Web: www.biotechsupportgroup.com