

Human Thrombin/Antithrombin Complex Antigen ELISA Kit

Catalog # HTATKT-COM

Strip well format. Reagents for up to 96 tests. Rev: July, 2017

REAGENTS PROVIDED

- •96-well antibody coated microtiter strip plate (removable wells 8x12) containing anti-human thrombin antibody, blocked and dried.
- •10X Wash buffer: 1 bottle of 50ml
- •Human TAT complex standard: 1 vial lyophilized standard
- •Anti-human antithrombin primary antibody: 1 vial lyophilized monoclonal antibody
- •Anti-mouse horseradish peroxidase secondary antibody: 1 vial concentrated HRP labeled antibody
- •TMB substrate solution: 1 bottle of 10ml solution

STORAGE AND STABILITY

Store all kit components at 4°C upon arrival. Return any unused microplate strips to the plate pouch with desiccant. Reconstituted standards and primary may be stored at -80°C for later use. Do not freeze-thaw the standard and primary antibody more than once. Store all other unused kit components at 4°C. This kit should not be used beyond the expiration date.

OTHER REAGENTS AND SUPPLIES REQUIRED

- •Microtiter plate shaker capable of 300 rpm uniform horizontally circular movement
- •Manifold dispenser/aspirator or automated microplate washer
- •Microplate reader capable of measuring absorbance at 450 nm
- Pipettes and Pipette tips
- •Deionized or distilled water
- Polypropylene tubes for dilution of standard
- Paper towels or laboratory wipes
- $\bullet 1N \ H_2SO_4 \ or \ 1N \ HCl$
- •Bovine Serum Albumin Fraction V (BSA)
- •Tris(hydroxymethyl)aminomethane (Tris)
- •Sodium Chloride (NaCl)

INTENDED USE

Human thrombin/antithrombin (TAT) complex antigen assay is intended for the quantitative determination of the covalent complex of thrombin and its inhibitor antithrombin III in human plasma and serum. **For research use only.**

BACKGROUND

Antithrombin III is a glycosylated plasma serine protease inhibitor that forms a stoichiometric complex with coagulation cascade enzymes [1]. Thrombin is a twochain vitamin K-dependent glycosylated serine protease that is activated from prothrombin in the coagulation cascade [2]. Antithrombin III inhibits thrombin with heparin enhanced kinetics and forms a 1:1 covalent complex [3].

ASSAY PRINCIPLE

Human TAT complex in samples will bind to the antihuman thrombin capture antibody coated on the microtiter plate. After appropriate washing steps, monoclonal anti-human antithrombin primary antibody binds to TAT complex captured on the plate. Excess antibody is washed away and bound monoclonal antibody is then reacted with the secondary antibody conjugated to horseradish peroxidase. TMB substrate is used for color development at 450nm. A standard calibration curve is prepared along with the samples to be measured using dilutions of TAT complex. Color development is proportional to the concentration of TAT complex in the samples. Free thrombin, prothrombin and antithrombin III will not be detected by this assay.

PRECAUTIONS

- •FOR LABORATORY RESEARCH USE ONLY. NOT FOR DIAGNOSTIC USE.
- Do not mix any reagents or components of this kit with any reagents or components of any other kit. This kit is designed to work properly as provided.
- •Always pour peroxidase substrate out of the bottle into a clean test tube. Do not pipette out of the bottle as contamination could result.
- •Keep plate covered except when adding reagents, washing, or reading.
- •DO NOT pipette reagents by mouth and avoid contact of reagents and specimens with skin.
- •DO NOT smoke, drink, or eat in areas where specimens or reagents are being handled.

PREPARATION OF REAGENTS

- •TBS buffer: 0.1M Tris, 0.15M NaCl, pH 7.4
- •Blocking buffer (BB): 3% BSA (w/v) in TBS
- •1X Wash buffer: Dilute 50ml of 10X wash buffer concentrate with 450ml of deionized water.

SAMPLE COLLECTION

Collect plasma using EDTA, heparin, or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000xg within 30 minutes of collection. Assay immediately or aliquot and store at \leq -20°C. Avoid repeated freeze-thaw cycles.

ASSAY PROCEDURE

Perform assay at room temperature. Vigorously shake plate (300rpm) at each step of the assay.

Preparation of Standard

Reconstitute standard by adding 1ml of blocking buffer directly to the vial and agitate gently to completely dissolve contents. This will result in a 1000ng/ml standard solution. Dilution table for preparation of human TAT complex standard:

TAT complex						
concentration	Dilutions					
(ng/ml)						
1000	From standard vial					
500	500μl (BB) + 500μl (1000ng/ml)					
250	500μl (BB) + 500μl (500ng/ml)					
100	600µl (BB) + 400µl (250ng/ml)					
50	500µl (BB) + 500µl (100ng/ml)					
25	500μl (BB) + 500μl (50ng/ml)					
10	600μl (BB) + 400μl (25ng/ml)					
5	500μl (BB) + 500μl (10ng/ml)					
2	600μl (BB) + 400μl (5ng/ml)					
1	500μl (BB) + 500μl (2ng/ml)					
0	500μl (BB) Zero point to					
0	determine background					

NOTE: DILUTIONS FOR THE STANDARD CURVE AND ZERO STANDARD MUST BE MADE AND APPLIED TO THE PLATE IMMEDIATELY.

Standard and Unknown Addition

Remove microtiter plate from bag and add 100µl TAT complex standards (in duplicate) and unknowns to wells. Carefully record position of standards and unknowns. Shake plate at 300rpm for 30 minutes. Wash wells three times with 300µl wash buffer. Remove excess wash by gently tapping plate on paper towel or kinwipe.

NOTE: The assay measures total human TAT complex in the 1-1000 ng/ml range. Samples giving human TAT complex levels above 1000ng/ml should be diluted in blocking buffer before use. Normal human plasma samples should be applied directly to the plate without dilution. A 1:500-1:1000 dilution for normal human serum is suggested for best results.

Primary Antibody Addition

Reconstitute primary antibody by adding 10ml of blocking buffer directly to the vial and agitate gently to completely dissolve contents. Add 100µl to all wells. Shake plate at 300rpm for 30 minutes. Wash wells three times with 300µl wash buffer. Remove excess wash by gently tapping plate on paper towel or kimwipe.

Secondary Antibody Addition

Briefly centrifuge vial before opening. Dilute 2µl of conjugated secondary antibody in 10ml of blocking buffer and add 100µl to all wells. Shake plate at 300rpm for 30 minutes. Wash wells three times with 300µl wash buffer. Remove excess wash by gently tapping plate on paper towel or kimwipe.

Substrate Incubation

Add 100 μ I TMB substrate to all wells and shake plate for 2-10 minutes. Substrate will change from colorless to different strengths of blue. Quench reaction by adding 50 μ l of 1N H₂SO₄ or HCl stop solution to all wells when samples are visually in the same range as the standards. Add stop solution to wells in the same order as substrate upon which color will change from blue to yellow. Mix thoroughly by gently shaking the plate.

Measurement

Set the absorbance at 450nm in a microtiter plate spectrophotometer. Measure the absorbance in all wells at 450nm. Subtract zero point from all standards and unknowns to determine corrected absorbance (A_{450}).

Calculation of Results

Plot A₄₅₀ against the amount of TAT complex in the standards. Fit a straight line through the linear points of the standard curve using a linear fit procedure if unknowns appear on the linear portion of the standard curve. Alternatively, create a standard curve by analyzing the data using a software program capable of generating a four parameter logistic (4PL) curve fit. The amount of TAT complex in the unknowns can be determined from this curve. If samples have been diluted, the calculated concentration must be multiplied by the dilution factor.

A typical standard curve (EXAMPLE ONLY):

Human TAT Complex Full Range 2.0 8 0.0 200 400 600 800 1000 0 TAT Complex (ng/mL) Human TAT Complex Linear Range 0.3 OD 450nm t= 6 minutes 0.2 0.1 0 0 10 20 30 40 50 TAT Complex (ng/mL)

EXPECTED VALUES

Concentration of TAT complex in normal human plasma was found to be 2.1 ng/ml and was increased in deep vein thrombosis (9.4 ng/ml) and disseminated intravascular coagulation (15.8 ng/ml) [4].

PERFORMANCE CHARACTERISTICS

Sensitivity: The minimum detectable dose (MDD) was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates (range OD_{450} : 0.099-0.117) and calculating the corresponding concentration. The MDD was 0.718 ng/ml.

Intra-assay Precision: These studies are currently in progress. Please contact us for more information.

Inter-assay Precision: These studies are currently in progress. Please contact us for more information.

Recovery: These studies are currently in progress. Please contact us for more information.

Linearity: These studies are currently in progress. Please contact us for more information.

Specificity: This assay recognizes natural human TAT complex. Natural human thrombin, prothrombin and antithrombin III were assayed for cross-reactivity. No significant cross-reactivity was observed.

Sample Values: Samples were evaluated for the presence of the antigen at varying dilutions.

Sample Type	Dilution	Mean (µg/mL)			
Citrate Plasma	Undiluted	0.01			
Serum	1:320	52			
Serum	1:640	53			

DISCLAIMER

This information is believed to be correct but does not claim to be all-inclusive and shall be used only as a guide. The supplier of this kit shall not be held liable for any damage resulting from handling of or contact with the above product.

REFERENCES

- 1. Rosenberg RD & Damus PS: J Biol Chem. 1973, 248:6490-6505.
- 2. Mann KG: Methods Enzymol. 1976, 45(Pt B):123-156.
- 3. Olson ST, et al.: Methods Enzymol. 1993, 222:525-559.
- 4. Hoek JA, et al.: Clin Chem. 1988, 34:2058-2062.

Example of ELISA Plate Layout 96 Well Plate: 22 Standard wells, 74 Sample wells

	1	2	3	4	5	6	7	8	9	10	11	12
Α	0	1	2	5	10	25	50	100	250	500	1000	
ŀ		ng/ml 1	ng/ml 2	ng/ml 5	ng/ml 10	ng/ml 25	ng/ml 50	ng/ml 100	ng/ml 250	ng/ml 500	ng/ml 1000	
В	0	ng/ml	ng/ml	ng/ml	ng/ml	ng/ml	ng/ml	ng/ml	ng/ml	ng/ml	ng/ml	
С								2				
D												
Ε							ッ		50			
F											S	
G												
н				5			V					

